


Anti-Rps27A Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62979

Item Number: CM0023833



Introduction

Rabbit polyclonal antibody against human RPS27A (ubiquitin-ribosomal protein eS31 fusion). Ideal for WB, IHC-P, IF/ICC, ELISA applications. Cross-reacts with human, mouse, rat. Validated in research on ubiquitin signaling and ribosome function.

[Learn More](#)

Parameter	Value
Product Name	RPS27A Rabbit pAb
Remarks/Alias	UBC; S27A; eS31; CEP80; UBA80; HEL112; UBCEP1; UBCEP80; RPS27A
Species	Human
Gene ID (Human)	6233
Gene ID	6233
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-156 of human RPS27A (NP_001170884.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62979
Protein Weight	18kDa
Shipping	Ice bag