


Anti-Rps28 Polyclonal Antibody For Wb, If/Icc, Elisa - P62857

Item Number: CM0022984



Introduction

High-quality rabbit anti-RPS28 polyclonal antibody validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat RPS28. Ideal for ribosome biogenesis and protein synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	RPS28 Rabbit pAb
Remarks/Alias	S28; eS28; DBA15; RPS28
Species	Human
GeneID (Human)	6234
GeneID	6234
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-69 of human RPS28 (NP_001022.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62857
Protein Weight	8kDa
Shipping	Ice bag