


Anti-Rps3 Polyclonal Antibody For Wb, If/Icc, Elisa - P23396

Item Number: CM0017578



Introduction

Rabbit polyclonal antibody against human RPS3, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Suitable for ribosomal protein detection and DNA repair studies.

[Learn More](#)

Parameter	Value
Product Name	RPS3 Rabbit pAb
Remarks/Alias	S3, uS3, RPS3
Species	Human
Gene ID (Human)	6188
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-243 of human RPS3 (NP_000996.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23396
Protein Weight	27kDa
Shipping	Ice bag