

Anti-Rps3 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P23396

Item Number: CM0010656



Introduction

High-quality rabbit monoclonal antibody against human RPS3 (SWISS: P23396). Validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for ribosomal protein research, DNA repair, and apoptosis studies.

[Learn More](#)

Parameter	Value
Product Name	RPS3 Rabbit mAb
Remarks/Alias	S3; uS3; RPS3
Species	Human
Gene ID	6188
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human RPS3 (P23396).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23396
Protein Weight	27kDa
Shipping	Ice bag