

Anti-Rps3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P23396

Item Number: CM0024016



Introduction

High-quality RPS3 Rabbit Polyclonal Antibody for WB, IHC-P, IF/ICC, ELISA. Reacts with Human, Mouse, Rat. Targets 40S ribosomal protein S3, involved in translation, DNA repair, and apoptosis. Ideal for cell biology research.

[Learn More](#)

Parameter	Value
Product Name	RPS3 Rabbit pAb
Remarks/Alias	S3; uS3; RPS3
Species	Human
GeneID (Human)	6188
GeneID	6188
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-243 of human RPS3 (NP_000996.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23396
Protein Weight	27kDa
Shipping	Ice bag