


Anti-Rps4X Rabbit Polyclonal Antibody For Wb And Elisa - P62701

Item Number: CM0026633



Introduction

Anti-RPS4X rabbit polyclonal antibody validated for WB and ELISA, reacting with human, mouse, and rat 40S ribosomal protein eS4 (30 kDa). Ideal for studying ribosome biogenesis and protein synthesis mechanisms.

[Learn More](#)

Parameter	Value
Product Name	RPS4X Rabbit pAb
Remarks/Alias	S4; eS4; CCG2; RPS4; SCAR; SCR10; DXS306; RPS4X
Species	Human
GeneID (human)	6191
GeneID	6191
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 66-263 of human RPS4X (NP_000998.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62701
Protein Weight	30kDa
Shipping	Ice bag