


Anti-Rps5 Rabbit Polyclonal Antibody For Wb And Elisa - P46782

Item Number: CM0030931



Introduction

Rabbit polyclonal antibody targeting human RPS5 (40S ribosomal protein S5, uS7), validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Ideal for ribosome biogenesis and protein synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	RPS5 Rabbit pAb
Remarks/Alias	S5; uS7; RPS5
Species	Human
GeneID (Human)	6193
GeneID	6193
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-204 of human RPS5 (NP_001000.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46782
Protein Weight	23 kDa
Shipping	Ice bag