


Anti-Rps6Ka1/Rps6Ka2/Rps6Ka3 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15418 / Q15349 / P51812

Item Number: CM0023263



Introduction

Rabbit polyclonal antibody against RPS6KA1, RPS6KA2, RPS6KA3, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for signal transduction and cancer research.

[Learn More](#)

Parameter	Value
Product Name	RPS6KA1/RPS6KA2/RPS6KA3 Rabbit pAb
Species	Human
GeneID (human)	6195/6196/6197
GeneID	6195/6196/6197
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 600-700 of human RPS6KA1/RPS6KA2/RPS6KA3 (NP_002944.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15418/Q15349/P51812
Shipping	Ice bag