

Anti-Rps6Ka2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q15349

Item Number: CM0022851

Introduction

Rabbit polyclonal antibody against RPS6KA2 (RSK3/p90RSK2), validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for probing MAPK/ERK signaling, cancer biology, and translational control.

[Learn More](#)

Parameter	Value
Product Name	RPS6KA2 Rabbit pAb
Remarks/Alias	RSK; HU-2; RSK3; p90RSK2; p90-RSK3; pp90RSK3; MAPKAPK1C; S6K-alpha; S6K-alpha2; RPS6KA2
Species	Human
Gene ID (Human)	6196
Gene ID	6196
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 632-733 of human RPS6KA2 (NP_066958.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15349
Protein Weight	83kDa
Shipping	Ice bag