

Anti-Rps6Ka4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 075676

Item Number: CM0014726

Introduction

Rabbit polyclonal antibody against ribosomal protein S6 kinase alpha-4 (MSK2, UniProt 075676). Validated for WB, IHC-P, ELISA; cross-reacts with human and mouse. 86 kDa. Ideal for inflammatory and stress signaling studies.

[Learn More](#)

Parameter	Value
Product Name	RPS6KA4 Rabbit pAb
Remarks/Alias	MSK2; RSK-B; S6K-alpha-4; RPS6KA4
Species	Human
Gene ID (Human)	8986
Gene ID	8986
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 673-772 of human RPS6KA4 (NP_003933.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	075676
Protein Weight	86kDa
Shipping	Ice bag