

Anti-Rps6Ka5 (Msk1) Polyclonal Antibody For Wb, Elisa - 075582

Item Number: CM0030978



Introduction

High-quality rabbit polyclonal antibody against human RPS6KA5 (MSK1), validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying kinase signaling in stress and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	RPS6KA5 Rabbit pAb
Remarks/Alias	MSK1; RLPK; MSPK1; RPS6KA5
Species	Human
GeneID (Human)	9252
GeneID	9252
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 700-800 of human RPS6KA5 (NP_004746.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075582
Protein Weight	90kDa
Shipping	Ice bag