


Anti-Rps6Ka5 (Msk1) Rabbit Monoclonal Antibody For Wb, Elisa - O75582

Item Number: CM0006370



Introduction

Rabbit monoclonal antibody targeting RPS6KA5 (MSK1), validated for WB and ELISA. Reacts with human, mouse, rat, and monkey. Ideal for studying mitogen- and stress-activated kinase signaling, transcription regulation, and histone modification pathways.

[Learn More](#)

Parameter	Value
Product Name	RPS6KA5 Rabbit mAb
Remarks/Alias	MSK1; RLPK; MSPK1; RPS6KA5
Species	Human
GeneID (Human)	9252
GeneID	9252
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 639-736 of human RPS6KA5 (NP_004746.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Monkey
SWISS	O75582
Protein Weight	62kDa/82kDa/90kDa
Shipping	Ice bag