

Anti-Rps6KL1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y6S9

Item Number: CM0027781

Introduction

Anti-RPS6KL1 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Recognizes human, mouse, and rat RPS6KL1, a ribosomal protein serine/threonine kinase. Ideal for ribosomal biology and signaling studies.

[Learn More](#)

Parameter	Specification
Product Name	RPS6KL1 Rabbit pAb
Remarks/Alias	RSKL2; RPS6KL1
Species	Human
Gene ID	83694
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human RPS6KL1 (NP_113652.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6S9
Protein Weight	60kDa
Shipping	Ice bag