

Anti-Rps8 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P62241

Item Number: CM0021253



Introduction

Rabbit polyclonal antibody against human RPS8 (40S ribosomal protein S8) for Western blot, IF/ICC, and ELISA. Validated in human, mouse, and rat samples. Predicted molecular weight 24 kDa.

[Learn More](#)

Parameter	Value
Product Name	RPS8 Rabbit pAb
Remarks/Alias	S8; eS8; RPS8
Species	Human
Gene ID (Human)	6202
Gene ID	6202
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-100 of human RPS8 (NP_001003.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62241
Protein Weight	24kDa
Shipping	Ice bag