


Anti-Rps9 Rabbit Polyclonal Antibody For Wb And Elisa - P46781

Item Number: CM0022643



Introduction

Rabbit polyclonal antibody targeting human RPS9 (40S ribosomal protein S9), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for ribosome biogenesis and translation research.

[Learn More](#)

Parameter	Value
Product Name	RPS9 Rabbit pAb
Remarks/Alias	S9; uS4; RPS9
Species	Human
GeneID (human)	6203
GeneID	6203
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-162 of human RPS9 (NP_001004.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46781
Protein Weight	23kDa
Shipping	Ice bag