


Anti-Rraga Polyclonal Antibody For Wb, Ihc-P, Elisa - Q7L523

Item Number: CM0012873



Introduction

Rabbit polyclonal RRAGA antibody validated for Western Blot (WB), Immunohistochemistry (IHC-P), and ELISA. Detects human, mouse, and rat RRAGA (FIP1/RAGA). Ideal for mTORC1 signaling and amino acid sensing research.

[Learn More](#)

Parameter	Value
Product Name	RRAGA Rabbit pAb
Remarks/Alias	FIP1; RAGA; FIP-1; RRAGA
Species	Human
GeneID (Human)	10670
GeneID	10670
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 213-313 of human RRAGA (NP_006561.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7L523
Protein Weight	37kDa
Shipping	Ice bag