


Anti-Rraga Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q7L523

Item Number: CM0030629



Introduction

Rabbit polyclonal antibody raised against human RRAGA (Ras-related GTP-binding protein A), validated for WB, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for mTORC1 signaling and amino acid sensing research.

[Learn More](#)

Parameter	Value
Product Name	RRAGA Rabbit pAb
Remarks/Alias	FIP1; RAGA; FIP-1; RRAGA
Species	Human
GeneID (Human)	10670
GeneID	10670
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 213-313 of human RRAGA (NP_006561.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7L523
Protein Weight	37kDa
Shipping	Ice bag