

Anti-Rras2 Rabbit Polyclonal Antibody For Wb, Elisa - P62070

Item Number: CM0029092



Introduction

RRAS2 Rabbit pAb is a polyclonal antibody validated for WB and ELISA, reactive against human, mouse, and rat. Ideal for studying Ras-related signaling and developmental processes.

[Learn More](#)

Parameter	Value
Product Name	RRAS2 Rabbit pAb
Remarks/Alias	NS12, TC21, RRAS2
Species	Human
Gene ID (Human)	22800
Gene ID	22800
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-204A of human RRAS2 (NP_036382.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62070
Protein Weight	23kDa
Shipping	Ice bag