

Anti-Rrbp1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9P2E9

Item Number: CM0013599

Introduction

Rabbit polyclonal antibody targeting human RRBP1, suitable for WB, IF/ICC, and ELISA applications. Cross-reacts with human and mouse. Ribosome-binding protein 1 mediates ribosome-ER membrane interaction.

[Learn More](#)

Parameter	Value
Product Name	RRBP1 Rabbit pAb
Remarks/Alias	RRp; hES; p180; ES130; ES/130; RRBP1
Species	Human
GeneID (Human)	6238
GeneID	6238
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human RRBP1 (NP_004578.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9P2E9
Protein Weight	152kDa
Shipping	Ice bag