

Anti-Rrm1 Monoclonal Antibody For Wb, Elisa - P23921

Item Number: CM0008661



Introduction

Recombinant rabbit monoclonal antibody targeting human RRM1, validated for Western blot and ELISA. Cross-reacts with mouse and rat. This 90kDa ribonucleotide reductase large subunit is essential for DNA synthesis, making the antibody ideal for oncology and cell cycle research.

[Learn More](#)

Parameter	Value
Product Name	RRM1 Rabbit mAb
Remarks/Alias	R1, RR1, RIR1, RRM1
Species	Human
GeneID (Human)	6240
GeneID	6240
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 691-792 of human RRM1 (P23921)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23921
Protein Weight	90kDa
Shipping	Ice bag