

Anti-Rrm1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P23921

Item Number: CM0012849



Introduction

Rabbit polyclonal antibody targeting human RRM1 (ribonucleotide reductase large subunit) for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with mouse and rat. Essential tool for studying DNA synthesis and related disorders.

[Learn More](#)

Parameter	Value
Product Name	RRM1 Rabbit pAb
Remarks/Alias	R1; RR1; RIR1; RRM1
Species	Human
GeneID (human)	6240
GeneID	6240
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 593-792 of human RRM1 (NP_001024.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23921
Protein Weight	90kDa
Shipping	Ice bag