

Anti-Rrm1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - P23921

Item Number: CM0015499



Introduction

Rabbit polyclonal anti-RRM1 antibody validated for WB, IF/ICC, IP, ELISA. Detects human and mouse ribonucleotide reductase large subunit (90kDa), essential for DNA synthesis and repair.

[Learn More](#)

Parameter	Value
Product Name	RRM1 Rabbit pAb
Remarks/Alias	R1; RR1; RIR1; RRM1
Species	Human
GeneID (human)	6240
GeneID	6240
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P23921
Protein Weight	90kDa
Shipping	Ice bag