


Anti-Rrm2 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P31350

Item Number: CM0029840



Introduction

Rabbit polyclonal antibody against human RRM2, suitable for WB, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with mouse and rat. Validated for research on nucleotide metabolism and DNA synthesis.

[Learn More](#)

Parameter	Value
Product Name	RRM2 Rabbit pAb
Remarks/Alias	R2; RR2; RR2M; C2orf48; RRM2
Species	Human
Gene ID (Human)	6241
Gene ID	6241
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-389 of human RRM2 (NP_001025.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31350
Protein Weight	45kDa
Shipping	Ice bag