

Anti-Rrm2 Rabbit Monoclonal Antibody For Wb And Elisa - P31350

Item Number: CM0011103



Introduction

Recombinant rabbit monoclonal antibody targeting human RRM2 (ribonucleotide reductase subunit M2), validated for WB and ELISA in human, mouse, and rat. Essential for DNA synthesis research.

[Learn More](#)

Parameter	Value
Product Name	RRM2 Rabbit mAb
Remarks/Alias	R2; RR2; RR2M; C2orf48; RRM2
Species	Human
GeneID (Human)	6241
GeneID	6241
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-117 of human RRM2 (P31350).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31350
Protein Weight	45kDa
Shipping	Ice bag