

Anti-Rrm2B Rabbit Polyclonal Antibody For Elisa And If/Icc - Q7Lg56

Item Number: CM0026866



Introduction

High-quality rabbit polyclonal antibody targeting human RRM2B (p53R2). Suitable for ELISA and IF/ICC. Ideal for DNA damage and repair research. Backed by CamelBio's comprehensive raw-material supply platform.

[Learn More](#)

Parameter	Value
Product Name	RRM2B Rabbit pAb
Remarks/Alias	P53R2; RCDFRD; MTDPS8A; MTDPS8B; RRM2B
Species	Human
GeneID (Human)	50484
GeneID	50484
Immunogen	synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RRM2B (NP_056528.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q7LG56
Protein Weight	41kDa
Shipping	Ice bag