


Anti-Rrp1 Rabbit Polyclonal Antibody For Wb And Elisa - P56182

Item Number: CM0020273



Introduction

High-quality rabbit polyclonal antibody targeting RRP1 (NNP-1, NOP52), suitable for Western blot and ELISA in human, mouse, and rat samples. Ideal for investigating 28S rRNA processing and nucleolar dynamics. SWISS: P56182.

[Learn More](#)

Specification	Details
Product Name	RRP1 Rabbit pAb
Remarks/Alias	NNP-1, NOP52, RRP1A, D21S2056E, RRP1
Species	Human
Gene ID	8568
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-253 of human RRP1 (NP_003674.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56182
Protein Weight	53kDa
Shipping	Ice bag