

## Anti-Rrp36 Polyclonal Antibody For Wb And Elisa - Q96Eu6

Item Number: CM0013804

### Introduction

Rabbit polyclonal antibody against human RRP36 (ribosomal RNA processing protein 36 homolog). Validated for Western blot and ELISA applications. Involved in pre-rRNA maturation leading to 18S rRNA.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | RRP36 Rabbit pAb                                                                                                 |
| Remarks/Alias    | C6orf153; dj20C7.4; RRP36                                                                                        |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 88745                                                                                                            |
| GeneID           | 88745                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-259 of human RRP36 (NP_149103.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human                                                                                                            |
| SWISS            | Q96EU6                                                                                                           |
| Protein Weight   | 30kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |