


Anti-Rrp8 Polyclonal Antibody For Wb, Elisa - 043159

Item Number: CM0023760



Introduction

Rabbit polyclonal antibody against RRP8 (Nucleomethylin) for Western blot and ELISA applications. Validated in human, mouse, and rat samples. Suitable for rDNA silencing and chromatin regulation studies.

[Learn More](#)

Parameter	Value
Product Name	RRP8 Rabbit pAb
Remarks/Alias	NML; KIAA0409; RRP8
Species	Human
GeneID (human)	23378
GeneID	23378
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-200 of human RRP8 (NP_056139.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43159
Protein Weight	51kDa
Shipping	Ice bag