

Anti-Rsad2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Wxg1

Item Number: CM0027807

Introduction

RSAD2 Rabbit pAb validated for WB, IF/ICC, ELISA with cross reactivity to human, mouse, rat. Interferon-inducible antiviral protein viperin; suitable for innate immunity and viral replication studies.

[Learn More](#)

Parameter	Value
Product Name	RSAD2 Rabbit pAb
Remarks/Alias	SAND; cig5; vig1; cig33; RSAD2
Species	Human
GeneID (Human)	91543
GeneID	91543
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 118-361 of human RSAD2 (NP_542388.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WXG1
Protein Weight	42kDa
Shipping	Ice bag