


Anti-Rsk2 / Rps6Ka3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P51812

Item Number: CM0003356



Introduction

Rabbit monoclonal RSK2/RPS6KA3 antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat RSK2 (84kDa). Essential for studying MAPK/ERK signaling, cell proliferation, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	RSK2 / RPS6KA3 Rabbit mAb
Remarks/Alias	CLS; RSK; HU-3; RSK2; MRX19; ISPK-1; XLID19; p90-RSK2; pp90RSK2; MAPKAPK1B; S6K-alpha3; RSK2 / RPS6KA3
Species	Human
GeneID (Human)	6197
GeneID	6197
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 600-700 of human RSK2 / RPS6KA3 (P51812).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51812
Protein Weight	84kDa
Shipping	Ice bag