


Anti-Rsl24D1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Uha3

Item Number: CM0021506



Introduction

Rabbit polyclonal antibody targeting human RSL24D1, a nucleolar protein involved in 60S ribosomal subunit biogenesis. Validated for WB and ELISA, reacts with human and mouse. Ideal for ribosome assembly studies.

[Learn More](#)

Parameter	Value
Product Name	RSL24D1 Rabbit pAb
Alias	L30; RLP24; RPL24; TVAS3; RPL24L; C15orf15; HRP-L30-iso; RSL24D1
Species	Human
GeneID	51187
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-80 of human RSL24D1 (NP_057388.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UHA3
Protein Weight	20kDa
Shipping	Ice bag