


Anti-Rspo1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q2Mka7

Item Number: CM0029086



Introduction

RSPO1 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, rat RSPO1. Ideal for investigating Wnt/ β -catenin pathway activation and stem cell niches. UniProt Q2MKA7.

[Learn More](#)

Parameter	Value
Product Name	RSPO1 Rabbit pAb
Remarks/Alias	RSPO; CRISTIN3; RSPO1
Species	Human
GeneID	284654
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 90-210 of human RSPO1 (NP_001033722.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q2MKA7
Protein Weight	29kDa
Shipping	Ice bag