

Anti-Rsrc1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96Iz7

Item Number: CM0026878

Introduction

Rabbit polyclonal antibody targeting human RSRC1 (SRrp53), validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for splicing and transcription research.

[Learn More](#)

Parameter	Value
Product Name	RSRC1 Rabbit pAb
Remarks/Alias	MRT70; BM-011; SFRS21; SRrp53; RSRC1
Species	Human
Gene ID (Human)	51319
Gene ID	51319
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 117-276 of human RSRC1 (NP_001258763.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96Iz7
Protein Weight	39kDa
Shipping	Ice bag