


Anti-Rsu1 Rabbit Polyclonal Antibody For Wb, Elisa - Q15404

Item Number: CM0017012



Introduction

Rabbit polyclonal antibody against human RSU1 (RSP-1), validated for WB and ELISA. Detects human, mouse, and rat RSU1. Suitable for Ras signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	RSU1 Rabbit pAb
Remarks/Alias	RSP-1; RSU1
Species	Human
GeneID (human)	6251
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-277 of human RSU1 (NP_036557.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15404
Protein Weight	32kDa
Shipping	Ice bag