

Anti-Rtcb Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y3I0

Item Number: CM0016057



Introduction

RTCB Rabbit Polyclonal Antibody validated for WB and ELISA; detects human, mouse, rat RTCB. Ideal for tRNA splicing and unfolded protein response research.

[Learn More](#)

Parameter	Value
Product Name	RTCB Rabbit pAb
Remarks/Alias	FAAP; HSPC117; C22orf28; DJ149A16.6; RTCB
Species	Human
GeneID (human)	51493
GeneID	51493
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y3I0
Protein Weight	55kDa
Shipping	Ice bag