


Anti-Rtf1 Rabbit Polyclonal Antibody For Wb, Elisa - Q92541

Item Number: CM0019282



Introduction

Rabbit polyclonal antibody against human RTF1 (GTL7, KIAA0252), a transcription elongation factor. Suitable for Western blot and ELISA. Detects human and mouse RTF1.

[Learn More](#)

Parameter	Value
Product Name	RTF1 Rabbit pAb
Remarks/Alias	GTL7, KIAA0252, RTF1
Species	Human
GeneID (Human)	23168
GeneID	23168
Immunogen	Synthetic peptide corresponding to amino acids 300-400 of human RTF1 (NP_055953.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q92541
Protein Weight	80kDa
Shipping	Ice bag