

Anti-Rtkin / Rhotekin Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bst9

Item Number: CM0027237

Introduction

Rabbit polyclonal antibody against human RTKN (Rhotekin), validated in WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Useful for studying Rho signaling, NF-kappa-B activation, and gastric tumorigenesis.

[Learn More](#)

Parameter	Value
Product Name	RTKN Rabbit pAb
Remarks/Alias	RTKN; rhotekin
Species	Human
GeneID (Human)	6242
GeneID	6242
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 254-513 of human RTKN (NP_001015056.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BST9
Protein Weight	63kDa
Shipping	Ice bag