


Anti-Rtn3 Rabbit Polyclonal Antibody For Wb And Elisa - O95197

Item Number: CM0016128



Introduction

Rabbit polyclonal antibody against human RTN3 (Reticulon-3) validated for WB and ELISA. Targets a protein involved in ER membrane trafficking, apoptosis, and viral response. Ideal for neuroscience and cell biology research.

[Learn More](#)

Parameter	Value
Product Name	RTN3 Rabbit pAb
Remarks/Alias	HAP; ASYIP; NSPL2; NSPLI; RTN3-A1; RTN3
Species	Human
GeneID (Human)	10313
GeneID	10313
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human RTN3 (NP_006045.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O95197
Protein Weight	113kDa
Shipping	Ice bag