


Anti-Rtn4 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9Nqc3

Item Number: CM0019302



Introduction

Rabbit polyclonal antibody against RTN4, validated for WB, IF-P, ELISA. Reacts with human, mouse, rat. Suitable for neurite outgrowth, ER tubule formation, and vascular remodeling studies. Ideal for diagnostic and research applications.

[Learn More](#)

Parameter	Value
Product Name	RTN4 Rabbit pAb
Remarks/Alias	ASY; NSP; NOGO; RTN-X; NSP-CL; RTN4-A; RTN4-C; RTN4-B1; RTN4-B2; NI220/250; Nbla00271; Nbla10545; RTN4
Species	Human
Gene ID (Human)	57142
Gene ID	57142
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human RTN4 (NP_065393.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NQC3
Protein Weight	130kDa
Shipping	Ice bag