

Anti-Rtn4Rl1 Rabbit Polyclonal Antibody For Wb, Elisa - Q86Un2

Item Number: CM0017668

Introduction

Rabbit polyclonal antibody against human RTN4RL1 (NgR3, NGRH2). Validated for WB and ELISA; cross-reacts with mouse and rat. Ideal for studying axon regeneration, neurite outgrowth inhibition, and neuroprotection pathways.

[Learn More](#)

Parameter	Value
Product Name	RTN4RL1 Rabbit pAb
Remarks/Alias	NgR3; NGRH2; RTN4RL1
Species	Human
Gene ID	146760 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-420 of human RTN4RL1 (NP_848663.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q86UN2
Protein Weight	49kDa
Shipping	Ice bag