

Anti-Rubicon Rabbit Polyclonal Antibody For Wb And Elisa - Q92622

Item Number: CM0024531



Introduction

High-quality rabbit polyclonal antibody against human Rubicon (RUBCN) for WB and ELISA applications. Targets an autophagy and immune regulatory protein, ideal for studying endosomal trafficking and pathogen defense mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Rubicon Rabbit pAb
Remarks/Alias	SCAR15; RUBICON; KIAA0226; Rubicon
Species	Human
GeneID	9711
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human Rubicon (NP_055502.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q92622
Protein Weight	109kDa
Shipping	Ice bag