

Anti-Rubicon Rabbit Polyclonal Antibody For Wb, Elisa - Q92622

Item Number: CM0015640



Introduction

Rabbit polyclonal antibody targeting human Rubicon (RUBCN) for Western blot and ELISA applications. Detects mouse and rat Rubicon. Recombinant immunogen corresponding to amino acids 1-200. Suitable for autophagy, endocytosis, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	Rubicon Rabbit pAb
Remarks/Alias	SCAR15; RUBICON; KIAA0226; Rubicon
Species	Human
GeneID (human)	9711
GeneID	9711
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human Rubicon (NP_055502.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q92622
Protein Weight	109kDa
Shipping	Ice bag