

Anti-Rufy1 Rabbit Polyclonal Antibody For Wb And Elisa - Q96T51

Item Number: CM0026681



Introduction

Rabbit polyclonal antibody against human RUFY1, suitable for Western blotting and ELISA. Detects mouse and rat RUFY1, a key regulator of endosomal cargo sorting and retrieval to the trans-Golgi network. Ideal for cellular trafficking research.

[Learn More](#)

Parameter	Value
Product Name	RUFY1 Rabbit pAb
Remarks/Alias	RABIP4; ZFYVE12; RUFY1
Species	Human
GeneID (human)	80230
GeneID	80230
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 439-708 of human RUFY1 (NP_079434.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q96T51
Protein Weight	80kDa
Shipping	Ice bag