

Anti-Rufy3 Rabbit Polyclonal Antibody For Wb And Elisa - Q7L099

Item Number: CM0015040

Introduction

Rabbit polyclonal antibody against RUFY3 (RIPX, SINGAR1, ZFYVE30) for Western blot (WB) and ELISA applications; suitable for human, mouse, and rat species; targets a 53 kDa protein.

[Learn More](#)

Parameter	Value
Product Name	RUFY3 Rabbit pAb
Remarks/Alias	RIPX; SINGAR1; ZFYVE30; RUFY3
Species	Human
GeneID (Human)	22902
GeneID	22902
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-65 of human RUFY3 (NP_055776.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q7L099
Protein Weight	53kDa
Shipping	Ice bag