

Anti-Runx1 Polyclonal Antibody For Wb, Elisa - Q01196

Item Number: CM0013044



Introduction

Rabbit polyclonal antibody targeting human RUNX1 transcription factor, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for hematopoiesis, leukemia, and immune research.

[Learn More](#)

Parameter	Value
Product Name	RUNX1 Rabbit pAb
Remarks/Alias	AML1; CBFA2; EVI-1; AMLCR1; PEBP2aB; CBF2alpha; AML1-EVI-1; PEBP2alpha; RUNX1
Species	Human
GeneID (Human)	861
GeneID	861
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human RUNX1 (NP_001001890.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01196
Protein Weight	49kDa
Shipping	Ice bag