


## Anti-Runx1T1 Polyclonal Antibody For Wb, Elisa - Q06455

Item Number: CM0024060



### Introduction

Rabbit polyclonal antibody against human RUNX1T1 (CBFA2T1/ETO/MTG8). Validated for WB and ELISA; cross-reacts with rat. Suitable for leukemia and transcriptional regulation studies.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | RUNX1T1 Rabbit pAb                                                                                                   |
| Remarks/Alias    | CDR; ETO; MTG8; AML1T1; ZMYND2; CBFA2T1; AML1-MTG8; RUNX1T1                                                          |
| Species          | Human                                                                                                                |
| GeneID (Human)   | 862                                                                                                                  |
| GeneID           | 862                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 250-403 of human RUNX1T1 (NP_783552.1) |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, ELISA                                                                                                            |
| Cross Reactivity | Rat                                                                                                                  |
| SWISS            | Q06455                                                                                                               |
| Protein Weight   | 68kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |