

Anti-Runx2 Rabbit Monoclonal Antibody For Wb, Ip, Elisa - Q13950

Item Number: CM0002818



Introduction

High-quality rabbit monoclonal antibody against human RUNX2, validated for Western blot, immunoprecipitation, and ELISA. Cross-reacts with mouse. Ideal for studying osteoblast differentiation and skeletal morphogenesis.

[Learn More](#)

Parameter	Value
Product Name	RUNX2 Rabbit mAb
Remarks/Alias	CCD; AML3; CCD1; CLCD; OSF2; CBFA1; OSF-2; PEA2aA; PEBP2aA; CBF-alpha-1; RUNX2
Species	Human
GeneID (Human)	860
GeneID	860
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13950
Protein Weight	57kDa
Shipping	Ice bag