

Anti-Runx3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q13761

Item Number: CM0028543



Introduction

Rabbit polyclonal antibody targeting human RUNX3 (AML2/CBFA3), validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Ideal for transcription regulation, cancer, and T cell development research.

[Learn More](#)

Parameter	Specification
Product Name	RUNX3 Rabbit pAb
Remarks/Alias	AML2; CBFA3; PEBP2aC; RUNX3
Species	Human
GeneID (Human)	864
GeneID	864
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 286-415 of human RUNX3 (NP_004341.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13761
Protein Weight	44kDa
Shipping	Ice bag