


Anti-Rusc1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9Bvn2

Item Number: CM0024219



Introduction

Anti-RUSC1 rabbit polyclonal antibody validated for WB, IF-P, ELISA. Reacts with human, mouse, rat RUSC1 (AP-4 complex accessory subunit, Nesca), involved in neuronal vesicular trafficking and differentiation.

[Learn More](#)

Parameter	Value
Product Name	RUSC1 Rabbit pAb
Remarks/Alias	NESCA; RUSC1
Species	Human
GeneID (Human)	23623
GeneID	23623
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human RUSC1 (NP_055143.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BVN2
Protein Weight	96kDa
Shipping	Ice bag