


Anti-Ruvbl1/Tip49A/Pontin Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y265

Item Number: CM0025968



Introduction

Rabbit polyclonal antibody targeting human RUVBL1/TIP49A/PONTIN (Q9Y265). Validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human and mouse. Ideal for chromatin remodeling and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	RUVBL1/TIP49A/PONTIN Rabbit pAb
Remarks/Alias	RVB1; TIH1; ECP54; TIP49; ECP-54; INO80H; NMP238; PONTIN; TIP49A; NMP 238; Pontin52; RUVBL1/TIP49A/PONTIN
Species	Human
GeneID (Human)	8607
GeneID	8607
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human RUVBL1/TIP49A/PONTIN (NP_003698.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y265
Protein Weight	50kDa
Shipping	Ice bag